

# RAPID Instrument Support

as of December 1, 2009

Prod Group	Model Names	PC Operating SW	Remote Image	Support PC	Comments/Limitations
<b>RAPID Supported</b>					
FIB	V600 Series	Windows 2000 Windows XP	Yes	Mandatory	
SDB	Helios Series	Windows XP	Yes	Mandatory	
SDB	Nova NanoLab Series	Windows 2000 Windows XP	Requires SW version XT V 2.0.5	Mandatory	
SDB	Strata Series	Windows 2000 Windows XP	Requires SW version XT V 2.0.5	Mandatory	
SDB	Quanta3D Series	Windows 2000 Windows XP	Yes	Mandatory	
SEM	Inspect Series	Windows 2000 Windows XP	Yes	Mandatory	
SEM	Nova NanoSEM Series	Windows 2000 Windows XP	Yes	Mandatory	
SEM	Quanta Series	Windows 2000 Windows XP	Yes	Mandatory	
TEM	Morgagni Series	Windows 2000 Windows XP	CCD cam required	Optional	Service cam not yet released for Morgagni
TEM	Tecnai Series	Windows 2000 Windows XP	Either CCD or Service cam required	Optional	For Tecnai10 and Tecnai Monochromator, V 2.1.5 is last release = Win2000
TEM	Tecnai Series	Windows NT	Either CCD or Service cam required	Mandatory	Windows XP recommended on support PC
TEM	Titan Series	Windows XP	Either CCD or Service cam required	Optional	
FPD	Certus3D Series	Windows 2000 Windows XP	Yes	Optional	
<b>RAPID Support Under Evaluation</b>					
FPD	Altura Series				
FPD	DB8xx Series				
FPD	Expida Series				
FPD	FIB 800 Series				
FPD	FIB 8xx Series				
FPD	XL 12xx Series				
FPD	XL 8xx Series				
SDB	DB235/237 (Strata)				
SDB	FIB 200 Series				
SEM	XL/Sirion Series				
<b>No RAPID Support</b>					
TEM	Vitrobot				
TEM	EMXXX Series				
TEM	CMXXX Series				
TEM	EM208 Series				
SEM	SEM 5xx Series				
SEM	Electro Scan Series				
FPD	XL50 Series				
CMD	All Series				
FPD	ATxxx Series				
FPD	FIB 610/611 Series				
FPD	FIB 820 & 825I Series				